



## IV Pro Adapter™

Mounting system for IV Pole / Accessory compatible with Safety Arm System™

### Characteristics

- Lightweight
- Perfectly fits clamp block shape
- Immediate access
- Match with Stretcher Safety Arm™ System
- Available both sides of the cot

### Application of use

- Ambulance
- EMS Vehicles
- Stretcher/cot
- Military Vehicles

### Advantages

- Minimum space for installation
- IV Pole keeps all its functions
- Resist to biological fluids
- Rust and sodium resistant
- 1 year full warranty

Any medical equipment that has been inadequately installed in ambulances, other emergency medical service (EMS) vehicles or crash carts, can potentially pose a serious risk or become dangerous in the event that the vehicle makes frequent or sudden stops, turns or rushed departures. This situation can result in damaging expensive medical equipment or even inflict injuries on a patient, ambulance personnel or clinic staff.



The **IV Pro Adapter™** is designed to allow the use of an IV pole holder on the same side as the stretcher **Safety Arm System™**. The IV Pro Adapter is lightweight and profiled to fit perfectly alongside the clamp block. With the addition of the adaptor, you now have the option to configure the stretcher to be equipped with one (1) or two (2) IV pole holders, located on opposite sides of the stretcher while maintaining the use of the Safety Arm System.

With this system, you gain the flexibility to place your IV pole holder to fit your operational needs. The adaptor is available to be installed on the left or right side of the stretcher (*looking from the foot end*). The installation process is very simple and takes only a few minutes. No additional hardware is required nor do you require to modify the IV Pole holder or the stretcher as you will use the existing hardware on the stretcher. Furthermore, the adaptor installation does not affect any of the safety aspects of the stretcher frame.

Once installed on the cot, the IV Pro Adapter allows you to use the IV Pole as in the original installation. It provides all the advantages of the original installation location, plus the capacity to also use the locking system of the clamp block of the SAS without any restrictions. It's the perfect accessory to maintain the location of the IV pole on the original stretcher configuration.

### About Us

At Technimount EMS, we understand that EMS and Critical Care Transport teams need solutions that they can rely on. We design and develop quality mounting systems which are safe and flexible, allowing you to concentrate on what matters most - patient care.

For use with the **Technimount System Safety Arm System™**

## **TECHNICAL SPECIFICATIONS**

<b>Product name:</b>	<b>IV Pro Adapter™</b>	
<b>Compatibility model:</b>	Technimount System – Safety Arm System Model 500 only Stryker cot – Serie 60 or 65	
<b>Environment of use:</b>	EMS / Hospital – Clinic – Biomed shop – Military hospital	
<b>Material:</b>	Plate: Aluminum high density w/ black anodize process	
<b>Color finish:</b>	Black velvet finish	
<b>Temperature of use:</b>	N/A	
<b>Impact Resistance<sup>1</sup>:</b>	<b>(30«G»)</b> Plus 10% safety factor <b>(On certified hook)</b>	
<b>Nbr of installation point:</b>	3 points on the original position on the cot	
<b>Dimensions (overall):</b>	15-3/16" Length X 3-5/8" Depth. X ¼" thickness (inch.)	
<b>Features:</b>	<ul style="list-style-type: none"> <li>- Available for right side or left side of the cot</li> <li>- Compatible with <b>Technimount cot system</b></li> <li>- Compatible with Stryker Cot</li> <li>- Free access to locking mechanism of the clamp block (SAS)</li> </ul>	
<b>Weight:</b>	Total (1 part) 1.10 lb.	
<b>Options:</b>	Available for the left side or right side (standard right side– Patient's leg) Color finish – contact sales for available colors	
<b>Locking mechanism:</b>	Original screw set from the Stryker Cot	
<b>Part number<sup>2</sup>:</b>	511-00-UN-RHS-KG	- Right side <b>with</b> Knee Gatch
	511-00-UN-RHS-WHKG	- Right side <b>without</b> Knee Gatch
	511-00-UN-LFS-KG	-Left Side <b>with</b> Knee Gatch
	511-00-UN-LFS-WHKG	-Left Side <b>without</b> Knee Gatch



### *Suggested Installation*

IV Pro Adapter, Technimount E.M.S are registered trademarks of Technologies CGC Inc

Stryker® and Stryker Cot® are trademarks of Stryker.

© 2021 Technimount E.M.S Holding Inc – All Rights Reserved - Printed in Canada.