


THE G3 CIRCULATORY KIT FEATURES A FOAM-LINED, FOLDABLE DESIGN THAT HELPS KEEP IV THERAPIES ORGANIZED AND AT YOUR GLOVED FINGERTIPS.



- Contents NOT included.
- Accommodates all IV bag sizes (1 Liter, 500 Milliliter, 250 Milliliter)
- Unique foam lined design keeps all Intravenous Emergency Medical Equipment inside protected from breakage or the elements
- Specifically designed to stow neatly in StatPacks Generation 3 Emergency Medical Services packs
- Hot or cold packs may be placed in to warm or cool fluid(s)
- Keeps IV Start Pack™ ready for use
- Front and inside pockets organize IV tubing and catheters together with the IV bag
- Transparent durable urethane window for quick IV bag solution identification
- Single-buckle closure securely holds contents in place and allows quick, one-handed access
- Top loop hangs the entire kit with the IV bag solution
- A quick stash pocket for trash or extra supplies
- Adjustable cam-lock buckle allows you to pressurize contents for rapid fluid infusion or high altitude fluid therapy
- Color-coded for easy identification of your EMS IV Kit



DRIP SETS, IV BAGS, CATHETERS, NEEDLES, AND START KITS ARE NEATLY ARRANGED FOR FAST ACCESS, AND DOUBLES AS A PRESSURE BAG FOR RAPID FLUID INFUSION.



THE G3 CIRCULATORY BAG HOLDS EVERYTHING YOU NEED FOR A QUICK IV SET-UP IN ONE CONVENIENT MODULE.



RED



BLACK

G3 Circulatory Kit

Contents NOT included

Accommodates all IV bag sizes (1 Liter, 500 Milliliter, 250 Milliliter)

Unique foam lined design keeps all Intravenous Emergency Medical Equipment inside protected from breakage or the elements

Specifically designed to stow neatly in StatPacks Generation 3 Emergency Medical Services packs

Hot or cold packs may be placed in to warm or cool fluid(s)

Keeps IV Start Pack™ ready for use

Front and inside pockets organize IV tubing and catheters together with the IV bag

Transparent durable urethane window for quick IV bag solution identification

Single-buckle closure securely holds contents in place and allows quick, one-handed access

Top loop hangs the entire kit with the IV bag solution

A quick stash pocket for trash or extra supplies

Adjustable cam-lock buckle allows you to pressurize contents for rapid fluid infusion or high altitude fluid therapy

Color-coded for easy identification of your EMS IV Kit

Size and Weight

Height: 10 inches

Width: 8 inches

Depth: 4 inches

Weight: 0.4 lbs

Available Colors

Red, Black

Capacity

Capacity: 210 in³

SKU

SKU: G36002RE [red]

SKU: G36002TK [black]