

**G3 VIAL STRANDS  
ACCOMMODATE MULTI-DOSE  
MEDS TIGHTLY IN PLACE WITH  
BRIGHT YELLOW PULL-TABS FOR  
EASY REMOVAL.**



- Contents NOT included
- Military spec Velcro™ hook on back attaches to any compartment with Velcro loop, which is a standard feature on the nearly all of our G3 packs, modules and cells. Each strand may hold 5 vials of various sizes
- No more cumbersome sleeves to access or difficulties identifying crucial medications



## G3 Vial Strand

Contents NOT included

Unique foam lined design keeps vials protected from breakage

Strong hook Velcro for secure attachment to any other loop Velcro

Specifically designed to stack neatly in StatPacks Generation 3 Emergency Medical Services packs

Bright yellow pulls for high visibility

3mm shock cord loops that stretch allowing vials to be neatly swapped out

Perfect for keeping random supplies together and organized

## Size and Weight

**Height:** 1.5 inches

**Width:** 5 inches

**Depth:** .25 inches

**Weight:** 0.12 lbs

## Available Colors

Red

## Capacity

Capacity: N/A

## SKU

SKU: G31098VS [Red]