

STATPACKS

THE G3 QUICKLOOK AED. THE ULTIMATE MOBILE ERGONOMIC FIRST AID STATION.

The G3 QUICKLOOK AED is built to accommodate most Automated External Defibrillators for visible and extremely prompt access. Perfect for use on rigs, personal vehicles and public areas such as airports, coliseums, shopping malls as an all-in-one first-responder AED backpack.



-  **ACLS** (Advanced Cardiovascular Life Support)
-  **AED** (Automated External Defibrillator)
-  **ALS** (Advanced Life Support)
-  **BVM** (Bag Valve Mask)
-  **PAL** (Pediatric Advanced Life Support)
-  **BLS** (Basic Life Support)
-  **O2** (Accommodates M6 or similar Oxygen Cylinder)



THE G3 QUICKLOOK AED EMPLOYS SIMPLE INTERNAL COMPARTMENTALIZATION THAT'S PERFECT FOR THE LAYPERSON AS WELL AS EXPERIENCED CLINICIANS.

Sussex County Emergency Medical Services has been using the Statpack Perfusion bags and medication modules for the past 6 years. During this period of time Statpacks has provided us with exceptional customer service, standing behind their products.

Bobby Schoonover, AA, NRP
Technical Services Division Manager - Sussex County, Delaware Emergency Medical Services



RED



BLUE



BLACK

QuickLook AED

- Contents NOT included
- A stow away "Rain Fly" is included for bad weather
- Unique design allows access to the main compartment from the backside of this Medical Backpack
- Single-buckle EMS AED pocket offers quick access to the Defibrillator
- Clear durable urethane window allows quick view of AED status
- Main compartment features three removable stick-in-place pouches that give visibility and organization to contents
- QuickZip access to main compartment through the back panel, maintaining access to contents and AED
- Tarpaulin bottom panel for protection from abrasion and the elements
- Foam lined shoulder straps that allow for comfortable transport of most emergency medical equipment
- Internal Velcro® lining strips to secure modules or other Velcro® compatible items to secure ambulance equipment
- Unique foam lined design keeps all Emergency Medical Equipment inside protected from breakage or the elements
- Side pockets for easy access to vital sign tools, bandages, and smaller ambulance equipment
- A quick stash pocket for trash and extra supplies located inside the main compartment
- ID window to label any way you want to distinguish between other EMS AED backpack's
- Main fabric is a non-toxic PVC-free tarpaulin. An environmentally friendly fabric featuring easy-clean, pathogen and water resistant properties

Size and Weight

- Height:** 19 inches
- Width:** 18 inches
- Depth:** 8 inches
- Weight:** 3.2 lbs

Available Colors

Red, Blue, Tactical Black

Capacity

Capacity: 1800 in3

SKU

- SKU: G35007RE [red]
- SKU: G35007BU [blue]
- SKU: G35007TK [tactical black]